



Hemodialysis Concentrates

Product Catalog



Driven to Deliver Life-Sustaining Dialysis Solutions™



Manufacturing & Delivering High-Quality Hemodialysis Concentrates to Dialysis Providers Worldwide

Rockwell Medical provides the hemodialysis community with high-quality¹ products, supported by a reliable delivery service, controlled by a robust Quality Management System, and regulated by the U.S. Food and Drug Administration (FDA). Rockwell Medical is ISO 13485 Certified and adheres to current Good Manufacturing Practices (cGMP) and Association for Advancement of Medical Instrumentation (AAMI) standards.

CitraPure®, Dri-Sate®, RenalPure®, and SteriLyte® are registered trademarks of Rockwell Medical, Inc.

¹Data on File. Available Upon Request. Source: Rockwell Medical, Inc.



Citric Acid Concentrate Powders and Liquids

CitraPure concentrate is citric-based and 100% acetate-free. CitraPure is formulated to be used in conjunction with RenalPure Powder Bicarbonate (CB-Series) or SteriLyte Liquid Bicarbonate (CL-Series) only in a 45X three-stream hemodialysis machine. The listed electrolyte values represent the final dialysate electrolyte yield.

45X CitraPure Liquid Acid Concentrate – Drum (55-gallon)

Formula (L4-Series)	Na+ (mEq/L)	K+ (mEq/L)	Ca++ (mEq/L)	Mg++ (mEq/L)	Chloride (mEq/L)	Citric Acid (mEq/L)	Dextrose (mg%)
L4-108-55	137.00	2.00	2.00	1.00	105.00	2.40	100.00
L4-110-55	137.00	3.00	2.00	1.00	106.00	2.40	100.00
L4-114-55	137.00	1.00	2.50	1.00	104.50	2.40	100.00
L4-116-55	137.00	2.00	2.50	1.00	105.50	2.40	100.00
L4-117-55	137.00	3.00	2.50	1.00	106.50	2.40	100.00
L4-124-55	137.00	2.00	3.00	1.00	106.00	2.40	100.00
L4-126-55	137.00	3.00	3.00	1.00	107.00	2.40	100.00

CitraPure, Dri-Sate, RenalPure, and SteriLyte are intended for use in acute and chronic hemodialysis and are to be used with a three-stream dialysis system. Custom formulas for 45X, 35X and 36.83X are available upon request.

Custom formulas are available upon request. Minimum volume requirements and premium pricing may apply.

CAUTION: Federal law restricts these devices to sale by or on the order of a physician.



45X CitraPure Liquid Acid Concentrate – Case (4 – 1 gallon (3.78 L) containers/case)

Formula (L4-Series)	Na+ (mEq/L)	K+ (mEq/L)	Ca++ (mEq/L)	Mg++ (mEq/L)	Chloride (mEq/L)	Citric Acid (mEq/L)	Dextrose (mg%)
L4-106-4	137.00	1.00	2.00	1.00	104.00	2.40	100.00
L4-108-4	137.00	2.00	2.00	1.00	105.00	2.40	100.00
L4-110-4	137.00	3.00	2.00	1.00	106.00	2.40	100.00
L4-114-4	137.00	1.00	2.50	1.00	104.50	2.40	100.00
L4-116-4	137.00	2.00	2.50	1.00	105.50	2.40	100.00
L4-117-4	137.00	3.00	2.50	1.00	106.50	2.40	100.00
L4-124-4	137.00	2.00	3.00	1.00	106.00	2.40	100.00
L4-126-4	137.00	3.00	3.00	1.00	107.00	2.40	100.00
L4-128-4	137.00	4.00	3.00	1.00	108.00	2.40	100.00



45X CitraPure Dry Acid Concentrate – Case (Each case produces 25 gallons of CitraPure Liquid)

Formula (D4-Series)	Na+ (mEq/L)	K+ (mEq/L)	Ca++ (mEq/L)	Mg++ (mEq/L)	Chloride (mEq/L)	Citric Acid (mEq/L)	Dextrose (mg%)
D4-103	137.00	2.00	2.25	1.00	105.25	2.40	100.00
D4-108	137.00	2.00	2.00	1.00	105.00	2.40	100.00
D4-110	137.00	3.00	2.00	1.00	106.00	2.40	100.00
D4-116	137.00	2.00	2.50	1.00	105.50	2.40	100.00
D4-117	137.00	3.00	2.50	1.00	106.50	2.40	100.00
D4-124	137.00	2.00	3.00	1.00	106.00	2.40	100.00
D4-126	137.00	3.00	3.00	1.00	107.00	2.40	100.00

CitraPure, Dri-Sate, RenalPure, and SteriLyte are intended for use in acute and chronic hemodialysis and are to be used with a three-stream dialysis system. Custom formulas for 45X, 35X and 36.83X are available upon request.

Custom formulas are available upon request. Minimum volume requirements and premium pricing may apply.

CAUTION: Federal law restricts these devices to sale by or on the order of a physician.



Acetic Acid Concentrate Powders and Liquids

RenalPure and Dri-Sate acid concentrates are formulated to be used in conjunction with RenalPure Powder Bicarbonate (CB-Series) or SteriLyte Liquid Bicarbonate (CL-Series) only in a 45X three-stream hemodialysis machine. The listed electrolyte values represent the final dialysate electrolyte yield.

45X RenalPure Liquid Acid Concentrate – Drum (55-gallon)

Formula (C-Series)	Na+ (mEq/L)	K+ (mEq/L)	Ca++ (mEq/L)	Mg++ (mEq/L)	Chloride (mEq/L)	Acetic Acid (mEq/L)	Dextrose (mg%)
C-108-55	137.00	2.00	2.00	1.00	105.00	4.00	100.00
C-110-55	137.00	3.00	2.00	1.00	106.00	4.00	100.00
C-114-55	137.00	1.00	2.50	1.00	104.50	4.00	100.00
C-116-55	137.00	2.00	2.50	1.00	105.50	4.00	100.00
C-117-55	137.00	3.00	2.50	1.00	106.50	4.00	100.00
C-124-55	137.00	2.00	3.00	1.00	106.00	4.00	100.00
C-126-55	137.00	3.00	3.00	1.00	107.00	4.00	100.00
C-130-55	137.00	2.00	2.25	1.00	105.25	4.00	100.00
C-131-55	137.00	3.00	2.25	1.00	106.25	4.00	100.00



45X RenalPure Liquid Acid Concentrate – Case (4 – 1 gallon (3.78 L) containers/case)

Formula (C-Series)	Na+ (mEq/L)	K+ (mEq/L)	Ca++ (mEq/L)	Mg++ (mEq/L)	Chloride (mEq/L)	Acetic Acid (mEq/L)	Dextrose (mg%)
C-102-4	137.00	0.00	0.00	1.00	101.00	4.00	100.00
C-106-4	137.00	1.00	2.00	1.00	104.00	4.00	100.00
C-107-4	137.00	2.00	1.50	1.00	104.50	4.00	100.00
C-108-4	137.00	2.00	2.00	1.00	105.00	4.00	100.00
C-110-4	137.00	3.00	2.00	1.00	106.00	4.00	100.00
C-114-4	137.00	1.00	2.50	1.00	104.50	4.00	100.00
C-116-4	137.00	2.00	2.50	1.00	105.50	4.00	100.00
C-117-4	137.00	3.00	2.50	1.00	106.50	4.00	100.00
C-120-4	137.00	1.00	3.00	1.00	105.00	4.00	100.00
C-124-4	137.00	2.00	3.00	1.00	106.00	4.00	100.00
C-126-4	137.00	3.00	3.00	1.00	107.00	4.00	100.00
C-130-4	137.00	2.00	2.25	1.00	105.25	4.00	100.00
C-131-4	137.00	3.00	2.25	1.00	106.25	4.00	100.00

CitraPure, Dri-Sate, RenalPure, and SteriLyte are intended for use in acute and chronic hemodialysis and are to be used with a three-stream dialysis system. Custom formulas for 45X, 35X and 36.83X are available upon request.

Custom formulas are available upon request. Minimum volume requirements and premium pricing may apply.

CAUTION: Federal law restricts these devices to sale by or on the order of a physician.



45X DriSate Dry Acid Concentrate – Case (Each case produces 25 gallons of RenalPure Liquid)

Formula (DC-Series)	Na+ (mEq/L)	K+ (mEq/L)	Ca++ (mEq/L)	Mg++ (mEq/L)	Chloride (mEq/L)	Acetic Acid (mEq/L)	Dextrose (mg%)
DC-108	137.00	2.00	2.00	1.00	105.00	4.00	100.00
DC-110	137.00	3.00	2.00	1.00	106.00	4.00	100.00
DC-116	137.00	2.00	2.50	1.00	105.50	4.00	100.00
DC-117	137.00	3.00	2.50	1.00	106.50	4.00	100.00
DC-124	137.00	2.00	3.00	1.00	106.00	2.40	100.00

Bicarbonate Concentrate Powders and Liquids

RenalPure Powder Bicarbonate and SteriLyte Liquid Bicarbonate are for use with CitraPure Dry or Liquid Acid Concentrate, Dri-Sate Dry Acid or RenalPure Liquid Acid Concentrate.



45X RenalPure Powder Bicarbonate

Formula (CB-Series)	Na+ (mEq/L)	Bicarbonate (mEq/L)	Product Description
CB-02-20C	37.00	37.00	1 bag yields 2.1 gallons. 20 bags/case
CB-15-4C	37.00	37.00	1 bag yields 15 gallons. 4 bags/case
CB-25-2C	37.00	37.00	1 bag yields 25 gallons. 2 bags/case



45X SteriLyte Liquid Bicarbonate

Formula (CL-Series)	Na+ (mEq/L)	Bicarbonate (mEq/L)	Product Description
CL-02-4	37.00	37.00	4 x 1 gallon jugs/case

CitraPure, Dri-Sate, RenalPure, and SteriLyte are intended for use in acute and chronic hemodialysis and are to be used with a three-stream dialysis system. Custom formulas for 45X, 35X and 36.83X are available upon request.

Custom formulas are available upon request. Minimum volume requirements and premium pricing may apply.

CAUTION: Federal law restricts these devices to sale by or on the order of a physician.



Concentrates Cleaner and Salt Pellets

Formula	Product Description
R-100-4	5% Acetic Acid Solution. 4 x 1 gallon/case.
PS-SALT25	25 lbs bag of salt. 1 Unit
PS-SALT50	50 lbs bag of salt. 1 Unit

Dry Acid Concentrate Mix System

1 Rapid Mixing

Mixes 100 gallons in approx. 3 minutes¹ – over 90% faster than traditional systems.

2 48" Angled Side Port


Ergonomic design for easy, mess-free powder loading, cleaning and maintenance.

3 Compact Footprint

Fits through standard doors. Ideal for space-conscious facilities.

4 Plug-and-Play Installation

Installs quickly. Simply connects to your clinic's existing power and RO water for consistent, high-quality performance.



5 One-Touch Operation

Simple interface. Just press the button – mixer stops automatically when ready to use.

6 No Onsite Technician Required

Self-installation and easy maintenance. No manufacturer service contracts required.

7 Made in the U.S.A.

Durable, corrosion-free stainless steel construction for long-term reliability. Fast turnaround from order to delivery.

Mixer Accessories

Item #	Description	Unit
PS-311	Hydrometer, Acid (C Series) Ranges 1.180 – 1.250	1 Unit
PS-24774-750	500 ml Hydrometer Cylinder	1 Unit
PS-GX05-9	Hytrex Replacement Cartridge Filter 5 micron Note: Replace when flow becomes weak	1 Unit
D4-Series or DC-Series	CitraPure Dry Acid Concentrate	25-gal mix
	Dri-Sate Dry Acid Concentrate	25-gal mix

Refer to Manufacturer's Instructions for Use for detailed technical specifications.

¹Data on File. Available Upon Request. Source: Rockwell Medical, Inc.

Committed to Improving the Lives of Dialysis Patients



Reliable & Responsive Customer Service is Our Priority

- Calculates minimum order requirements to help avoid extra fees
- Verifies all codes to help prevent ordering errors
- Creates transportation schedule to ensure consistent and timely product delivery
- Selects most cost effective delivery option for urgent orders (common carrier, drop ship, etc.)



Dedicated Transportation Services Ensure Delivery of Your Vital Orders

- Date and time for delivery is provided
- Tracks all orders and deliveries
- Deliver product to customer's desired storage location[†]
- Fleet is equipped with satellite tracking equipment
- Pump-off services for drum users[†]
- Drivers remove and recycle empty drums[†]

[†]Additional fees may apply

Driven to Deliver Life-Sustaining Dialysis Solutions™



Sales

sales@rockwellmed.com

Customer Care

(800) 449-3353

custserv@rockwellmed.com

Shipping and Deliveries

(248) 487-9085

shipping@rockwellmed.com

www.rockwellmed.com